

K8 CLINICAL CONNECTIONS & SETTINGS GUIDE

K8 allows pattern recognition command of Steeper Espire Pro powered elbow flexion/extension with optional pronation/supination of powered wrists and pattern recognition command of multifunction grips for various* hands. (*consult the Coapt Compatibility Guide)

PART 1: PHYSICAL CONNECTIONS


PRE-ASSEMBLED
The Coapt controller (white PCB) will be pre-attached to the Espire unit. Do NOT attempt to remove it.

PRE-ASSEMBLED CONNECTION
The Coapt **Complete Calibrate** button is attached to the Coapt controller with its flex cable passing through the Espire Strain Relief Disc. Contact Coapt support for assistance detaching /attaching this cable from the elbow.

PRE-ASSEMBLED CONNECTION
Ensure the 2-pin, white Coapt aux connection is attached

EMG CONNECTION
The Coapt EMG Interface Cable (not shown here) for Espire has a special 90° exit. Ensure it is connected and secured once the Espire Strain Relief Disc is in place.

STEEPER ESPIRE PRO & MULTIFUNCTION HAND



DISTAL CONNECTIONS

HANDS: TASKA, TASKA CX, COVVI, ZEUS, PSYONIC, (older) BEBIONIC

WRIST: Fillauer Standard Wrist Rotator (S103043, S103054, S103055) with 6-Band Coaxial Plug (S103059)

COLOR	OUTPUT
Orange	Hand Open
Green	Hand Close
Blue	Rotation
Red	Power
White	Hand Open
White	Hand Close

May need to be connected in reverse polarity when using COVVI

HAND: I-LIMB

WRIST: Fillauer Standard Wrist Rotator (S103043, S103054, S103055) or Ottobock Electric Wrist Rotator (R029)

COLOR	OUTPUT
Orange	Hand Open
Green	Hand Close
Blue	Rotation
Red	Power
White	Hand Open
White	Hand Close
Orange	Open Digital
Green	Close Digital

HAND: VINCENTEVOLUTION4

WRIST: Fillauer Standard Wrist Rotator (S103043, S103054, S103055) or Ottobock Electric Wrist Rotator (R029)

COLOR	OUTPUT
Orange	Hand Open
Green	Hand Close
Blue	Rotation
Red	Power
White	Hand Open
White	Hand Close
Orange	Open Digital
Green	Close Digital

This Digital Provided by Coapt

PART 2: SOFTWARE SETTINGS

ELBOW

STEEPER ESPIRE PRO

In the Espire iOS user application, most screens will indicate the direct Coapt connection and not require specific settings.

WRISTS

OTTOBOCK ELECTRIC WRIST ROTATOR (NO MYOROTRONIC), MOTION CONTROL STANDARD WRIST ROTATOR, MOTION CONTROL STANDARD WRIST ROTATOR & 6-BAND COAXIAL PLUG (NO PROCONTROL)

No settings required

HANDS

I-LIMB QUANTUM* (*WITH SERIAL NUMBERS M9618 OR NEWER)

Ensure hand settings are the factory defaults. For proper operation of Gesture Control via Coapt, ensure 'HO' is selected as the i-mo trigger in the i-Limb app

PSYONIC ABILITY HAND

Ensure the "Control System" setting in the PSYONIC ABILITY HAND app is set to "Coapt".

Ensure all other hand settings in the PSYONIC ABILITY HAND are set to factory defaults.

AETHER ZEUS HAND

Input Option: EMG, Input Sites: Dual Direct, Speed Control Strategy: Proportional, Grip Switching Mode: Open-Open

HANDS (CONT.)

OTTOBOCK PAST GENERATION BEBIONIC

Contact Coapt if needed

VINCENTEVOLUTION4

SENSOR SETUP:

Enable Sensor 1 (Open) - On, Enable Sensor 2 (Close) - On, Enable Sensor 3 - Off, Enable Sensor 4 - Off

Red Signal (Open): Threshold - 5, Gain - 30, Activation: 15

Blue signal (Close): Threshold - 5, Gain - 30, Activation: 15

Enable first signal wins strategy - Off

Time settings: Hold open - 3.0s max, Timeout - 7.0s min

Standard trigger settings: no trigger

MULTICHANNEL TRIGGER CONTROL:

Enable Multichannel Trigger Control - on.

Standard Trigger - No Trigger

Multichannel Hold Signal 1: lateral grasp with hold open

Multichannel Hold Signal 2: pinch grasp with hold open/close

Multichannel Hold Signal 3: index grasp with hold close

Multichannel Hold Signal 4: disabled

Period Peak For Multichannel Peaks: 0.2s

Multichannel Peak 1: cylinder grasp with close peak

Multichannel Peak 2: cupholder grasp with open peak

Multichannel Peak 3: power grasp with open peak/close peak

Multichannel Peak 4: disabled

COVVI HAND

Ensure all "Triggers" in the Grip Switch Trigger settings in the COVVI Go app are set to "Trigger Off"

Ensure all other hand settings in the COVVI Go app are set to the factory defaults



TASKA/TASKA CX

Settings will automatically be set upon each power-on. Use QR code for more information in the Taska|Coapt Clinician Guide.

