


# 13

## CLINICAL CONNECTIONS & SETTINGS GUIDE

I3 controls pronation/supination of powered wrists and direct pattern recognition command of multifunction grips for i-Limb Quantum\* or VINCENTevolution3. (\*with serial numbers M9618 or newer)

### PART 1: PHYSICAL CONNECTIONS



### PART 2: SOFTWARE SETTINGS

#### WRISTS

##### OTTOBOCK ELECTRIC WRIST ROTATOR (NO MYOROTRONIC)

No settings required

##### MOTION CONTROL STANDARD WRIST ROTATOR (NO PROCONTROL)

No settings required

#### HANDS

##### I-LIMB QUANTUM\*

Ensure hand settings are the factory defaults.

For proper operation of Gesture Control via Coapt, ensure 'CC' is selected as the i-mo trigger in the i-Limb app

##### VINCENTEVOLUTION3

Ensure hand settings are the factory defaults