


KO

# CLINICAL CONNECTIONS & SETTINGS GUIDE

KO allows pattern recognition command of Ottobock DynamicArm Standard elbow flexion/extension with optional pronation/supination of a powered wrist and open/close operation of prosthetic hands. Pattern recognition command of multifunction grips is not included.

## PART 1: PHYSICAL CONNECTIONS




## PART 2: SOFTWARE SETTINGS

### ELBOW

#### OTTOBOCK DYNAMICARM STANDARD

Ensure all settings in the Ottobock Data Station 'Customizing' screen are set exactly as shown.



### HANDS (CONT.)

#### OTTOBOCK SENSORHAND, MYOHAND VARIPLUS, SYSTEM ELECTRIC GREIFER, OR SYSTEM ELECTRIC HAND

Ensure 'control mode'/'program' 1

#### MOTION CONTROL PROHAND, PROETD/ETD2

##### HAND SET UP

- Input Type: EMG
- Hand Filtering: Quick
- Channel: Dual
- Control: Differential
- Hand Direction: A

##### Thresholds/Outputs

- A: 15
- B: 15

##### Input Gains

- A: 5
- B: 5

### WRISTS

#### OTTOBOCK ELECTRIC WRIST ROTATOR (NO MYOROTRONIC)

No settings required

### HANDS

**\*\*\*MANY/MOST HANDS SHOULD BE SET TO FACTORY DEFAULT FOR OPTIMAL OPERATION.**

#### PSYONIC ABILITY HAND

Ensure the "Control System" setting in the PSYONIC ABILITY HAND app is set to "Direct Control"

Ensure all other hand settings in the PSYONIC ABILITY HAND app are set to the factory defaults

#### OTTOBOCK BEBIONIC

Control Strategy: Dual Site

Open / Close Strategy: Electrodes: A,B

Control Options: Electrodes

Control Response: Proportional

In ELECTRODE SCREEN:

- Channel 1 "Blue" threshold: 10%
- Channel 1 "Red" threshold: 90%
- Channel 2 "Blue" threshold: 10%
- Channel 2 "Red" threshold: 90%

Remember to press "Send to Hand" after making any changes