

K2

CLINICAL CONNECTIONS & SETTINGS GUIDE

K2 allows pattern recognition command of LTI Boston Digital Arm elbow flexion/extension with optional pronation/supination of powered wrists and open/close operation of prosthetic hands. Pattern recognition command of multifunction grips is not included.

PART 1: PHYSICAL CONNECTIONS



CONNECTION

Coapt cable with circular elbow-mating connector attaches firmly to oriented connections on the elbow

BOSTON DIGITAL ELBOW
with Open/Close Terminal Device

NOTE:

When using an i-Limb Quantum terminal device, connections for **Open** and **Close** (at the distal end of the forearm) may need to be reversed to ensure proper operation.

PART 2: SOFTWARE SETTINGS

ELBOW

LTI BOSTON DIGITAL ARM

The Boston Digital Arm should be supplied directly from LTI and pre-setup with correct internal programming and cabling for Coapt use. Contact LTI support with any questions.

WRISTS

OTTOBOCK ELECTRIC WRIST ROTATOR (NO MYOROTRONIC)

No settings required

MOTION CONTROL STANDARD WRIST ROTATOR (NO PROCONTROL)

No settings required

HANDS

I-LIMB, TASKA, VINCENTEVOLUTION3, PSYONIC ABILITY HAND, COVVI HAND

Ensure hand settings are the factory defaults.

OTTOBOCK SENSORHAND, MYOHAND VARIPLUS, SYSTEM ELECTRIC GREIFER, OR SYSTEM ELECTRIC HAND

Ensure 'control mode'/'program' 1

HANDS (CONT.)

OTTOBOCK BEBIONIC

Control Strategy: Dual Site

Open / Close Strategy: Electrodes: A,B

Control Options: Electrodes

Control Response: Proportional

In ELECTRODE SCREEN:

- Channel 1 "Blue" threshold: 10%
- Channel 1 "Red" threshold: 90%
- Channel 2 "Blue" threshold: 10%
- Channel 2 "Red" threshold: 90%

Remember to press "Send to Hand" after making any changes

MOTION CONTROL PROHAND, PROETD/ETD2

HAND SET UP

- Input Type: EMG
- Hand Filtering: Quick
- Channel: Dual
- Conatrol: Differential
- Hand Direction: A

Thresholds/Outputs

- A: 15
- B: 15

Input Gains

- A: 5
- B: 5