


# EMG INTERFACE CABLE OPTIONS

## EMGE: STANDARD, RING TERMINALS AND FULL SET OF CONTROLSEAL™ DOME ELECTRODES


17 ring terminations for use with Coapt ControlSeal™ Dome Electrodes

\*Replacement or Additional Dome Electrode set available as p/n ELSB



## EMGS: EVALUATION-STYLE SNAPS

17 snap terminations for use with set of liner-snap electrode contacts (sold by third party) or with disposable, self-adhesive snap contacts (p/n ELEV).



## EMGC: EVALUATION-STYLE CUFF

Wrap-around, Velcro cuff for evaluation use with transradial and small circumference transhumeral limbs. Not suitable for definitive use.



## EMGO: EMG INTERFACE CABLE, OSSEOINTEGRATION-STYLE CUFF

Elastic cuff for definitive Osseointegration use only.

\*Sizing information may be requested when ordering.



**OTHER EMG INTERFACE CABLE OPTIONS MAY BE AVAILABLE UPON SPECIAL REQUEST.**