


H4

CLINICAL CONNECTIONS & SETTINGS GUIDE

H4 controls direct pattern recognition command of multifunction grips for the Nexus Hand by COVVI where no powered wrist or elbow is present.

PART 1: PHYSICAL CONNECTIONS



PART 2: SOFTWARE SETTINGS

HANDS

COVVI NEXUS HAND

Ensure all "Triggers" in the Grip Switch Trigger settings in the COVVI Go app are set to "Trigger Off"

Ensure all other hand settings in the COVVI Go app are set to the factory defaults