


18

CLINICAL CONNECTIONS & SETTINGS GUIDE

18 controls pronation/supination of Motion Control "6-band" powered wrist and direct pattern recognition command of multifunction grips for the Nexus Hand by COVVI.

PART 1: PHYSICAL CONNECTIONS



PART 2: SOFTWARE SETTINGS

WRISTS

MOTION CONTROL STANDARD WRIST ROTATOR & 6-BAND COAXIAL PLUG

No settings required

HANDS

NEXUS

Ensure hand settings are the factory defaults